10 Implants sur mesure: Custom-made implants

Avantages

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Imaging protocol for the manufacturing of custom hip prosthesis

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CONVENTIONNAL RADIOGRAPH

Teleradiography of 2 legs support
Installation of the subject:

- Parallel feet and extended knee
- The pelvis and last lumbar vertebrae must appear last on the radiograph
- In the event that there is an inequality of the leg length, a second telemetry must be performed with pelvic balance by a wedge placed under the shorter limb.
- The thickness of this block should be marked on the plate and allow full extension of the two knees
- If there is a doubt about the scale we recommend a grid layout (dimension known cases) or association with a graduated marker placed in the plane of the femurs.

Radiograph of the hip to operate

It is necessary to place a metal marker of known size in the plane of the femur

- 2 impacts: front and profile

The ideal is to receive emails or We-Transfer of digital radios

The images must contain the greatest possible diaphysis (minimum 5 cm longer than the length of the stem to manufacture)

1) facial image

- The femur must face (strict)

In some cases (stiffness, this can lead to place the subject with a scale marker in the plane of the femur)

- Be included on the image:

- the lower part of the sacra-iliac joint
- the outer edge of the greater trochanter
- the symphysis pubis

2) Profile image

- Strict profile of the femur
- The best defining is the proximal epiphysis and metaphysis

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1) Installation of the subject:

- Strict Immobility (wedging in rotation)
- Lower member parallel to the edge of the table
- patella at the zenith
- Sagittal wedge for horizontal proximal femur

2) Sequence programming on Front topogramme:

(All cuts must be numbered)

- *The first cut through the center of the head (C1)*
- *Definition first sequence (neck-shaft)*
- Minimum: 5 mm
- height: C1 at the base of the small trochanter (CT)
- *Definition of the second sequence (shaft)*
- Minimum: 10 mm
- height: C1 to CD, diaphyseal cut located at 180 mm from the first cut
- In the case where the femur would already be instrumented (presence of a femoral stem, plate osteosynthesis or others ...) The last cut (CD) as to be located to 50 mm in minimum below the equipment in place
- *Definition third sequence (helitortion)*
- 1 section through the condylar tubers, overlays two collar cuts and the posterior condylar plan pour calculation anteversion.

3) Acquisition protocol

- Thickness of cut: 1 mm
- Acquisition incremental volumic for three-dimensional reconstruction

Information required concerning the patient:


- The name, first name, date of birth, weight, height, side.

ADDITIONAL INFORMATION


In the event that there is a length difference of the lower limbs, it is desirable, to erase these measurement errors, to know the functional length inequality measured in this lying subject between the anterior superior iliac spine and a bone relief ankle.

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Zone de coupe scanners
Cutting Area Scanners

Angle cervico-daiphysaire (face)
Cervico-daiphyseal angle (face)

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




Fabrication de la prothese
Manufacture of the prosthesis

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


Hip arthroplasty in Pediatric Patient with Cerebral Palsy

Pediatric custom made total hip prosthesis

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


Pelvis modelling

The scanner cuts are imported into Osirix. Image processing steps allow us to isolate the patient's pelvis and create a mesh file, which can then be converted into a volume in SolidWorks.


This model is then cut according to the surgeon's instructions, depending on the surgical procedure.

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

Design of the prosthesis

The design starts from the cut of the pelvis given by the surgeon. The upper face of the prosthesis has to fit the shape of the pelvis cut. Finally, two anchoring stems are added to the upper face; an impacting hole and a cavity for the hip prosthesis are dug in the lower face.

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

Design of the instruments

The guide's shape is obtained by following the same procedure as for the prosthesis, which means by extruding matter from the cut of the pelvis. The guide is then drilled twice, at locations corresponding to the anchoring stems of the prosthesis. These holes are designed for the insertion of to drill bit guides which will allow the surgeon to drill the pelvis at the location of setting up of the stems.


Assembly

Once all the parts of the prosthesis are designed, they can be assembled and inserted into the patient's radiography. The prosthesis drawing had to be approved by the surgeon before being sent to the production line. At this step, 3D printings of the bone and the prosthesis can be performed.

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Custom made nail fixation and cementoplasty of the iliac bone

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
Référence
50.020.70

The supra-acetabular part of the iliac bone transmits weighting forces from the hip to the sacro-iliac joint and spine. Osteolytic metastases limited to the supra-acetabular area can be treated by cementoplasty which increases compression bone strength. However, when bone destruction extends to the arcuate line of the iliac bone and toward the sacroiliac joint, where shearing forces are exerted, cementoplasty is not sufficient to restore bone strength. In such situations, a nail introduced obliquely through the posterior inferior iliac spine toward the supra-acetabular region, along the thickest part of the iliac bone, is able to substitute the sustain forces from the hip to the sacro-iliac joint.

The nail is stabilized by fixing it by two screws placed posteriorly into the sacral ala through the sacroiliac joint and one screw placed into the anterior part of the nail into the supra-acetabular region. Finally, polymethylmetacrylate cement is injected through the nail for further bone fixation.

Indication : Extensive tumoral destruction of the supra-acetabular iliac bone including the arcuate line with preservation of the hip joint space and sacral ala.

Contraindication: The use of iliac nail instrumentation is contraindicated in patients with open wounds, infection, pregnancy, and osteolytic tumors destroying the hip joint space and sacral ala.


Référence	Longueur / Length
36.760.20	20
36.760.25	25
36.760.30	30
36.760.35	35
36.760.40	40
36.760.45	45
36.760.50	50
36.760.55	55
36.760.60	60
36.760.65	65
36.760.70	70
36.760.75	75
36.760.80	80
36.760.85	85

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The design and manufacturing process are carried out in five main steps:

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Before starting the design, the 3-D bone model was transferred to parametric software. The osteolytic bone lesion of a patient are analyzed. Followed by the definition of the dimensions and shape of the nail.

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3-D bone reconstruction of the iliac and sacrum bones.

The CT images are properly processed to generate the bone contours and create the patient-specific 3-D bone model

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
A drawing of the system and a 3-D printed prototype are created, and the design is approved by the surgeon.

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An initial validation of a system motion simulation, (implant and instrumentation), is performed (simulation of the nail implantation and screw locking mechanism). A Finite Element Method simulation is done in order to ensure the mechanical properties of the system.

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Imaging protocol for manufacture of custom knee prosthesis

Teleradiography of the 2 lower limbs under load Position of the subject: - feet together parallel and knees extended - the pelvis and the last lumbar vertebrae must appear on the image - in the case of an inequality in length, a second telemetry must be carried out with balance of the pelvis by a wedge placed under the shortest limb. - The thickness of this block must be noted on the photo and allow the complete extension of both knees - If there is any doubt about the scale, we recommend a grid photo (dimension of the known boxes) or the association with a graduated mark placed in the plane of the femurs. X-ray of the knee to be operated on.

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2 incidences: face and profile The ideal being to receive digital radios and the scanner via Wettransfert. The images must contain as much of the diaphysis as possible (at least 5 cm more than the length of the tibial and femoral stems to be manufactured).

Installation of the subject:

- strict immobility (rotational locking) - lower limb parallel to the edge of the table - ball joint at the zenith

- Programming the sequences on the front topogram: (all cuts must be numbered) - the first cut passes through the center of the head (C1) - definition of first sequence (cervico-diaphyseal) • increment: 5 mm • height: C1 at the base of the lesser trochanter (CT) - second sequence definition (diaphysis) • increment: 10 mm • height: C1 to CD, diaphyseal cut located 180 mm from the first cut • in the case where the bone is already instrumented (presence of a femoral stem, an osteosynthesis plate or others, etc.) the last cut (CD) must be located at least 50 mm below the material in place - definition third sequence (helitortion) 1 cut passing through the tubercles condylars, superpositions of two neck sections and the posterior condylar plane for calculation of anteversion.

Acquisition protocol - thickness of sections: 1 mm - Incremental volume acquisition for three-dimensional reconstruction Information to be provided concerning the patient: - name, first name, date of birth, weight, height, side

In the event that there is an inequality in the length of the lower limbs, it is desirable to know the inequality in functional length measured in this subject lying between the anterior superior iliac spine and a bony relief of the ankle.

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Radiographic images of a male patient who underwent a total knee replacement for an osteosarcoma of the left distal femur.

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Pre-operative radiograph showing a lytic lesion of the distal metaphysis and epiphysis of the left femur.

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Pre-operative MRI confirming that the epiphysis is involved with the tumor.

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Post-operative radiograph after tumor resection and total knee replacement with cemented megaprosthesis.

Bibliographie/ **Bibliography**

• Knee function after limb salvage surgery for malignant bone tumor: comparison of megaprosthesis and distal femur allograft with epiphysis sparing.

Sébastien Pesenti, Emilie Peltier, Vincent Pomero, Guillaume Authier, Lionel Roscigni, Elke Viehweger, Jean-Luc Jouve
International Orthopaedics (SICOT) DOI 10.1007/s00264-017-3608-x, August 2017.

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Stage 1: Preoperative planning and radiographic evaluation

Before the implant placement, the surgeon realize a preoperative planning and a radiographic evaluation which allow to preset the device type and the approximately size which appears the most suitable to the patient's morphology. During this crucial stage, the surgeon can be detected potential anatomy abnormalities in the patient and therefore anticipated potential complications which may compromise the successful completion of the surgical intervention.

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


Stage 2: Trial patellar implant implantation

Once the knee is exposed and resected, place the trial patellar implant on the exposed patella to choose the most suitable size. Then, this trial device is used to improve the device's placement and orientation. Remove the trial implant and implant the definitive patellar implant.

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
Planning a left tibial Growth prosthesis

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
Growing Knee prosthesis for malignant tumor under arthroscopy

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The progress of chemotherapy justifies surgical techniques which improve the survival quality of those stricken with osteosarcomas and Ewing's sarcomas. If the reconstruction techniques are reliable for adults, for children they run into difficulties because of the problem of future disparity of the limbs's length, particularly the reach of the femur's inferior extremity.

This localisation represents more than forty percent of osteosarcomas and close to twenty percent of Ewing's sarcomas. It is because of these cases that we suggest the growing prosthesis, as soon as the foreseeable disparity, due to the suppression of the growth cartilage of the inferior femur, is equal to or greater than 3 cm.

This prosthesis may be used right away after the tumoral resection, or placed secondarily when the disparity starts to interfere in functionality.

Description

The prosthesis is composed of two parts connected by a hinge. Only a hinged joint can avoid the luxation that lengthening would produce.

The tibial piece is not cemented, so as not to sterilize the growth cartilage of the tibia. The smooth stem is small enough to reduce growth cartilage trauma. When setting, perforation of the growth cartilage should be perfectly centered to limit deviations in the Valgus Varus. The tibial plate is equipped with two small anti-rotation flanges, meant to be included in the trench cut out of the articular cartilage and sub-chondral bone above the conjugation cartilage. The tibial plate rests upon the tibial articular cartilage.

The femoral part created to-size is equipped with a telescopic system that you can operate from the interior of the knee under arthroscopy, thus limiting as much as possible the cutaneous opening and the risk of infection (principal complication of other growing prostheses).

The particularly effective anti-twisting system avoids the releasing of metallic particles when pressure is exerted on the prosthesis.

The tail of the femoral prosthesis, according to need, may be cemented or screwed. It should be solidly set.

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Used for positioning the prosthesis' axle. The axle has a screw-foot so that one may hold it with a screwdriver

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Ancillary material

Trial axle

Before cementing the femoral piece, it is prudent to verify the length of the limb at the end of the intervention, and for this, use the trial axle. This prosthesis allows one to make up for an over-cutting error peroperatively, by lengthening the prosthesis. Once you have verified the correct length of the limb, you can cement the femoral implant, and assemble the tibial and femoral fragments with the rings and axle. To avoid the axle coming out of place, it is imperative that you place an anti-extrusion clip with clip applying forceps and verify that the clip is well in place, by turning it within its housing. Closing may then follow; do not hesitate to review the transfer of the medial twin when the musculo-cutaneous covering of the implant seems uncertain.

Lengthening technique

This is performed, if possible, after the end of chemotherapy, when the ailing person has recovered normal immunary defense. The arthroscopy is done according to the usual anaesthesia technique on a knee at 40° flexion. The arthroscopy indicates the intercondylar indentation of the prosthesis in which is located the head of the growing screw. Inserted through a 3 mm incision, the screw-driver must be pushed all the way into the screw-head's recess before being turned.

In order to begin the lengthening, turn the screw-driver counter-clockwise, if it is a question of a left or right prosthesis. When the lengthening becomes too long, the stretched quadriceps muscle tends to lengthen the knee and diminish the flexion. One must not try to lengthen on a knee bent at less than 30°, because the results will be painful and distinguished by a persistent flexion. Respecting this rule also avoids vasculo-nervous complications observed with other models.

When growth has stopped, it is advised to change the growing prosthesis, whose size is no longer adapted to that of the adolescent, to an adult reconstruction model prosthesis.

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


Planification de genou unicompartmentale Unicompartmental Knee Planning

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


The geometry of the knee is modelled in SolidWorks using the patient's data. The bone is cut according to the surgeon's instructions.

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The prosthesis is designed in order to reproduce the patient's knee geometry ; stems are added to the femoral and to the tibial parts to allow the attachment to the bone ; the femoral part is hollowed out so as to lighten the device.

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


Design of the instruments

In order to place the prosthesis correctly, the surgeon will be guided thanks to various instruments. The first of them is a cutting guide : it is applied to the femur as it reproduces its shape, is fixed thanks to $\varnothing 2$ pins, and allows the surgeon to cut the femur along its landmarks. Once this is done, a drilling guide can be applied to the cut area : it comprises two holes, which have the same orientation as the two femoral stems of the prosthesis. This guarantees a perfect fit of the prosthesis in the patient's knee.

The same procedure is then applied to the tibial part.

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