
11

Divers • Various

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Benefit of an unguial prosthesis and setting technique. Dr S. Guero

1. In the Emergency Room


Avulsion of the unguial tablet is accompanied by an impairing of the nail bed with two unsatisfactory developments: hyper-budding and keratinization.

Hyper-budding


The nail bed being very vascularized, botryomycomas are frequent. They are sometimes resistant to all medical treatment (corticosteroid, nitration) and, in some cases, surgical excision is necessary. Replacing the tablet or placing an unguial prosthesis mechanically prevents hyper-budding.

Keratinization

Henceforth, we know that the nail bed is not a «sterile matrix» (Fig. 1). According to the individual, keratin production varies between one and ten percent of the thickness of the unguial tablet. In the case of an unguial avulsion, after the scarring-over phase, the nail bed reacts to air exposure by starting to produce keratin. In the best case, the newly-formed nail is supple and butts upon the hardened nail bed. The nail slows its progress from proximal to distal and forms a cushion proximal to the nail bed. This cushion will only then thicken very slowly. The unguial tablet will arrive at the hyponychium in six months, instead of three or four, and will be very unattractive. Six more months must pass before it becomes refined.

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2. "Cold" surgery for unguial dystrophies

Consultation requests for unguial dystrophies are very frequent, since the esthetics of the hands and feet are important to patients. Of traumatic, infectious, or congenital origin, unguial dystrophies may be improved, even healed, provided that the regrowth of the new nail is guided in order to avoid physiological mishaps such as the ones mentioned above. Because this removal is painful, prescribe numerous bandages and sometimes cause unguial dystrophies at the given site. The remedy of an unguial prosthesis, whose make-up is close to the nail morphology, simplifies considerably the method of this surgery. The unguial prostheses are thus indispensable in repairing surgery:

- in urgent cases where the extracted nail is no longer available,
- "cold", in order to avoid the removal of a healthy nail.

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Bibliographie/ Bibliography :

- **Stéphane Guero** : Benefits of an Ungual Prosthesis in Traumatology and Reconstructive Surgery of the Nail
Techniques in Hand & Upper Extremity Surgery, Volume 18, Number 1, March 2014

Technical note


The ungual prosthesis is slightly larger than normal in order to enable resizing and adaptation to all nail forms. Its concavity reproduces the form of a normal nail. The fundamental quality is that its underlying side is striated longitudinally like a natural nail. The nail must be attached laterally by two separate points, or with a U-point (Fig. 2a and Fig. 2b).

On the other hand, it is not advised to place a point in the dorsal skin at the risk of arriving in the matrix and leaving a permanent longitudinal ridge (Fig. 2c). This setting technique advocated by certain authors must be done away with.

Moreover, the prosthesis will only be inserted under the skin during "cold" surgery, in reconstructing the groove of the matrix.

In urgent case, we suggest inserting a gap between the dorsal skin and the sub-cutaneous part of the nail (soft nail), which decrease the risks of periungual panaritium. Even if a prosthesis delivered sterile is better tolerated than an avulsed nail, it is imperative to insure proper drainage of the hæmatoma.


Échelle 1
Scale 1

Implant stérile
Sterile implant

Figure 2c

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La cupule d'interposition inox

Stainless steel Interposition cupula

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Cupula prostheses are interposition prostheses between the metatarsal head and the re-cut phalanx, set to the capsule with catgut or 000 nylon stitches.

These prostheses do not have an intra-phalangeal fixing eye-bolt, thus not causing bone reactions, such as osteolysis. Being thin (8/10 millimeters), they not cause intolerance or inflammatory reactions.

Acting as a meniscus, they prevent bone ankylosis and insure the continuation of articularity amplitude.

They are especially appropriate in surgical therapy of arthroses like Hallus Rigidus, Freiberg Koehler, and in iatrogenic and post-traumatic arthroses.

Interphalangeal use of the Regnault Cupula

In the case of a forty-years-old man:
painful stiffening, unaligned proximal inter-phalanx of the second toe after traumatic dislocation, reduced in emergency, unimmobilized.

Operation at 4 months : arthroplasty of the proximal inter-phalanx with stringy interposition.

On the fourth day, direct shock to the second toe, the results are very painful with toe distension, progressive stiffening, appearance of capsular ossification, pain caused by the least mobilization of the second toe.

Six months after the operation : painful stiffness of the proximal inter-phalanx.

The patient wants to regain mobility ; to avoid the impaction of P2 on P1, the interposition of a Regnault cupula re-cut to the dimensions proximal inter-phalanx is planned, after the recut of P1 and P2, and to exerece the capsular ossifications.

Operation realised in september 1990.


Six months after the second operation, the second toe has recovered its normal volume, exhibits no pain, flexion is satisfactory, extension is normal. Only a few vascular troubles persist. The cupula will be removed only if it poses a problem.

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

Mèche non canulée Ø 2 mm, Lg totale 50 mm, Lg taillée 18 mm ref 61.330.18
Non- cannulated Ø2 mm drill, total lg 50 mm, Cuted lg 18 mm

- Lame-plaque en titane pour ostéotomie de la 1^{re} phalange du gros orteil (raccourcissement ou varisation)
Blade-plate in titanium for osteotomy of the first phalanx of the hallux (shortening or varition)

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Indication

The titanium button is used for trans-osseous sutures of the cuff to support the node when the bone is soft. It may be used for all trans-osseous sutures by U-threading in which the bone is soft (fractures of the trochanter, tuberosity, etc.).

Thread to use

Any nylon thread, but it must be sufficiently resistant so as not to be quickly cut in two by the button.

Why titanium ?


In order to permit modern post-operative imaging via a scanner or M.R.I.

Noteworthy: The width of the button for the shoulder is very small, so there is no conflict with the trochanter and no protrusion.

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Contrôles: J 45, 6 mois, 1 an

For a long time, the treatment of osteoid osteoma involved a "block" osseous excision that was damaging and disproportionate with respect to the size of the lesion.

The difficulties in locating the nidus during surgery and the desire to limit the bone resection as much as possible contributed to the development of the percutaneous excision technique with scanographic location and guidance.

Advantages of PEBORD

Precision of the location

- Very limited nature of the excision (not very invasive procedure)
- Little risk of post-operative fracture
- Short-term hospitalization (48 to 72 hours)
- Rapid recovery of ability to use legs (with canes)

Resumption of normal activity is possible 2 months after the treatment

Percutaneous resection technique: PEBORD

- General anesthesia
- Collaboration of surgeon and radiologist
- Duration of the procedure: 45 to 90 minutes (locating and resection)
- Hospitalization: 1 to 2 days
- Immediate partial support

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1 Accurate locating of the nidus

The examination is done under local or general anesthesia. The lesion is located by making meeting tomographic cuts in order to define the ballistic strategy. A cut centered on the nidus is made to estimate its diameter and depth. In order to facilitate the drilling, the member is positioned so that the angle of attack is as direct as possible (as vertical as possible).

2 Percutaneous resection

A two centimeter-long cutaneous incision is made at the chosen point of entry marked on the skin. The following steps are then carried out:


- 1) The trocar with its mandrel is placed facing the lesion (image 1)
 - 2) The mandrel is replaced by a pin, which must go into the nidus. The correct position of the pin is verified by several tomographic cuts. (image 2)
 - 3) Once the pin has been properly positioned facing the nidus, the trocar is introduced up to the surface of the bone. A cover is slipped onto the trocar, which is then removed. (image 3 and 4)
 - 4) A 6 mm diameter drill bit is introduced on the pin. The drilling is done with a motorized unit. The same action is then done with a 9 mm diameter drill bit. (image 5)
 - 5) After the drilling, the zone of the nidus is approached by the trephine, which extracts the bone sample containing it. (images 6 and 7)
- Then milling is necessary to destroy the adjacent region to the nidus to prevent a possible relapse.
- 6) After the procedure, a scanographic check is done to verify complete excision and the absence of complications. If the excision is incomplete, a new drilling is done.

Few severe complications have been reported and the technique is less expensive than surgical resection.

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
Extraction kit of an established prosthesis

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
The trocar allows the drilling and the injection of osseous marrow.

The trocar allows a percutaneous use with an entrance point in the trochanteral area, an aiming of the femoral head under amplifier of brightness, and an intraosseous penetration with the hammer.

When the trocar is placed, one can leave in place the sleeve and withdraw the chuck. This sleeve then allows the injection of osseous marrow or another liquid substance using a syringe of 10 mL.

Then : reinstallation of the chuck and ablation.

- This device allows the sampling and the injection of marrow in same operational time.

Instrument for making punctures composed of :


- a metallic socket
- a mandrel : metallic rod that has a triangular tip with sharp edges at its end that goes 10 mm beyond the socket.

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