7

Prothèses • Prostheses

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Concept : The use of a long anatomic prosthesis, connected by the bank bone, seems to be the best solution for restoring the upper extremity of the femur when it is destroyed for oncological, orthopedic or traumatological reasons.

Indications : Iterative repair of cemented hip prostheses. The destruction of the bone mass caused by successive and iterative interventions on hip prostheses can require either the implanting of a reconstruction prosthesis in the remaining bone coupling, or the complete reconstruction of a larger or smaller part of the upper extremity of the femur by a reconstruction prosthesis connected by the allotransplant.

This method can increase the bone mass and decrease the mass of usable cement.

In the event of use of a prosthesis connected by an allotransplant, we recommend not cementing the part of the prosthesis in contact with the allotransplant, but simply the distal part fixed in the receiving bone.

Traumatology : Complex fracture of the upper extremity of the femur on arthritic, necrotic or destroyed head. The sealing is made on the distal part of the prosthesis. The proximal portion will be coupled to the autologous bone fragments.

Oncology : Wide surgical resection of the upper extremity of the femur and the adjacent muscles for bone tumor.

The use of an allotransplant connecting the prosthesis allows for rapid refixing of the peripheral muscles and reinsertion of the gluteus medius by transosseous points.

Advantages : Allows the patient to resume activity after a few days, allows for rapid refixing of the peripheral muscles, eliminates the peri-prosthetic chamber and decreases the leverage and increases the bone mass and puts the patient back in conditions of basic surgery.

Technical specificities : Anatomical implant with forged titanium alloy support TA6V ELI inducing a 10° anteversion with 3 mm lateralization for a cervicodiaphysal angle of 135°.

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Fond de cotyle de reconstruction Reconstruction acetabular support

Référence gauche / left	Référence droit / right	Ø en mm
50.084.46	50.083.46	46
50.084.48	50.083.48	48
50.084.50	50.083.50	50
50.084.52	50.083.52	52
50.084.54	50.083.54	54
50.084.56	50.083.56	56
50.084.58	50.083.58	58
50.084.60	50.083.60	60

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Surgical technique for trial acetabular cup and femoral head

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Stage 1: Preoperative planning and radiographic evaluation

Before the implant placement, the surgeon realize a preoperative planning and a radiographic evaluation which allow to preset the device type and the approximately size which appears the most suitable to the patient's morphology. During this crucial stage, the surgeon can be detected potential anatomy abnormalities in the patient and therefore anticipated potential complications which may compromise the successful completion of the surgical intervention.

Stage 2: Acetabulum preparation

The surgical approach for the cup's implantation can be varied according to the surgeon's habits. For any approach chosen, it is essential to have a perfect exposition of the acetabulum. In the first time, excise the acetabular labrum and if necessary, next realize a residuals osteophytes removal to obtain a good visualization of the acetabulum and to facilitate the drilling.

Stage 3: Cup's drilling

The cup's drilling need to be very progressive using increasing size of drills (2 to 2mm). The objective is to create an anatomic shape's cavity which ensure a good anchoring of the cup in a well vascularized osseous support.

Stage 4: Trial cup implantation

Assemble the trial cup which had the same size that the last used drill, on the replaceable handle for trial cup. Position the trial cup in the acetabulum and verify the size and adaptation to the cavity. If necessary, a complementary drilling can be realized.

Stage 5: Definitive cup implantation

Definitive cup's size is the same that the last used drill. Ensure that the patient is in a good position. Place the definitive cup on the acetabular cup impactor and implant it on the acetabulum cavity.

Stage 6: Trial femoral head implantation

Place the trial femoral head on the femoral stem and ensure that the selecting size, fits perfectly with the chosen cup. Remove the trial head and implant the definitive femoral head.

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Standard modular shoulder prosthesis

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